


Materials of construction

TangenX® SIUS® Gamma TFF Device components

1. TangenX® SIUS® Gamma Cassette
2. Cassette gaskets
3. Tri-clover clamps
4. Pinch clamps
5. Braided tubing
6. CPC AseptiQuik® connectors
7. Hose barb fittings
8. Tri-clover gaskets
9. Filter plate insert



Component	Material of construction
Cassette gaskets	Ethylene propylene diene rubber
Tri-clover clamps	Nylon glass filled
Pinch clamps	Nylon glass filled
Braided tubing	Platinum cured silicone
CPC AseptiQuik® connectors	Polycarbonate
Hose barb fittings	Polypropylene
Tri-clover gaskets	Platinum cured silicone
Filter plate insert	Polypropylene

TangenX® SIUS® Gamma Cassette materials of construction

- Modified polyether sulfone membrane
- Polypropylene / Polyethylene substrate
- Polypropylene screen
- High density polyethylene channel spacer
- Platinum cured silicone pressure sensitive adhesive
- Polyurethane adhesive

